

2 × AceTaq® Master Mix

P411/P412



Version 9.1

Vazyme biotech co., Ltd.

Introduction

This product includes AceTaq® DNA Polymerase, dNTPs, and an optimized buffer system that allows for amplification by simply adding primers and templates, reducing pipetting and improving throughput and reproducibility. A protective agent added to the amplification system allows the 2 × AceTaq® Master Mix to maintain stable activity after repeated freeze-thaw cycles. This product is available in a version containing electrophoresis buffer and dye. It can be directly electrophoresed after the reaction. The 3' end of the PCR product, A, can be directly cloned into the T vector and used in the ClonExpress® Cloning Kit (C112/C113/C114).

Components

Components	P411-01	P411-02	P411-03
2 × AceTaq Master Mix	1 ml	5 × 1 ml	15 × 1 ml

Components	P412-01	P412-02	P412-03
2 × AceTaq Master Mix (Dye Plus)	1 ml	5 × 1 ml	15 × 1 ml

Storage

Stored at -20°C.

Protocol

Reaction system

ddH ₂ O	To 50 µl
2 × AceTaq Master Mix	25 µl
Primer1 (10 µM)	2 µl
Primer2 (10 µM)	2 µl
Template DNA*	x µl

*The optimal reaction concentration of different templates is different. The following table shows the recommended template usage for 50 µl reaction system.

Human genomic DNA	1 - 500 ng
E.coli genomic DNA	1 - 100 ng
λDNA	0.1 - 1 ng
Plasmid DNA	0.1 - 1 ng

Reaction procedure

95°C	5 min (Pre-denaturation) ^a	
95°C	30 sec	} 30-35 cycles
55°C ^b	30 sec	
72°C	60 sec/kb	
72°C	7 min (Complete Extension)	

a. The pre-denaturation time takes at least 5 minutes. If the amplification is not ideal, the pre-denaturation time of 95°C can be extended appropriately, up to 10 min.

b. The annealing temperature needs to be adjusted according to the T_m value of the primer, and is generally set to be lower than the primer T_m value of 3 - 5°C.



Vazyme Biotech Co., Ltd.
www.vazyme.com

Order: global@vazyme.com Support: support@vazyme.com
For research use only, not for use in diagnostic procedures.

Notes

Primer Designing

1. Choose C or G as the last base of the 3' end of the primer;
2. Avoid continuous mismatch at the last 8 bases of the 3' end of the primer;
3. Avoid hairpin structure at the 3' end of the primer;
4. Tm of the primers should be between 55°C - 65°C;
5. 5' adding sequence should not be included when calculate Tm of the primers;
6. GC content of the primers should be between 40% - 60%;
7. Tm and GC content of forward and reverse primers should be as similar as possible.



ISO 9001: 2015