

# AceTaq DNA Polymerase

Catalog# P401

Version 9.1



Vazyme biotech co., ltd.

## Introduction

AceTaq DNA Polymerase is a chemically modified Taq DNA Polymerase that is completely blocked at room temperature and is released only after heating at 95°C. Non-specific amplification and primer dimerization can be prevented during sample preparation and reaction ramps. Compared with the antibody-based hot-start Taq, the polymerase activity of AceTaq is blocked more stringently and completely. It takes only 5 min to activate AceTaq. AceTaq is compatible with most existing PCR protocols. Combined with an optimized buffer system, AceTaq minimizes non-specific amplification and primer dimers, ensuring extremely high sensitivity and specificity, which make it ideal for amplifying low-copy genes from complex templates. The 3'-end of the PCR products contain A, which can be directly cloned to a T-vector and is suitable for ClonExpress and Topo Cloning Kits (#C112 / #C113 / #C115 / #C601).

## Components

Components	P401-d1	P401-d2	P401-d3
10 × AceTaq Buffer (Mg <sup>2+</sup> Plus)	1 ml	4 × 1 ml	
dNTP Mix (10 mM each)	200 µl	800 µl	3 × P401-d2
AceTaq DNA Polymerase (5 U / µl)	50 µl	200 µl	

## Storage

Store at -30°C ~ -15°C. Transport at -20°C ~ 0°C.

## Unit Definition

One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble material in 30 min at 74°C, using activated salmon sperm DNA as template/primer.

## Protocol

### Reaction system

ddH <sub>2</sub> O	to 50 µl
10 × AceTaq Buffer (Mg <sup>2+</sup> plus)	5 µl
dNTP Mix (10 mM each)	1 µl
Template DNA*	x µl
Primer 1 (10 µM)	2 µl
Primer 2 (10 µM)	2 µl
AceTaq DNA Polymerase (5 U / µl) **	0.5 µl

\* The optimal concentration is different for different templates. The following table shows the recommended template usage for a 50-µl reaction system:

Human genomic DNA	1-500 ng
<i>E. coli</i> genomic DNA	1-100 ng
λ DNA	0.1-1 ng
Plasmid DNA	0.1-1 ng

\*\* The amount of AceTaq can be adjusted between 0.25 µl-1 µl. Increasing the amount of enzyme may increase the yield of amplification, but it may also reduce the specificity.



Vazyme Biotech Co., Ltd  
www.vazyme.com

Order: [global@vazyme.com](mailto:global@vazyme.com)

Support: [global@vazyme.com](mailto:global@vazyme.com)

**For research use only, not for use in diagnostic procedures.**

## Reaction Program

95°C	5 min *(pre-denaturation)	
95°C	30 sec	} 30 -35 cycles
55°C **	30 sec	
72°C	60 sec / kb	
72°C	7 min (complete extension)	

\* The pre-denaturation takes at least 5 min. If the amplification is not ideal, extend the pre-denaturation up to 10 min.

\*\* The annealing temperature needs to be adjusted according to the T<sub>m</sub> value of the primer (normally set as be 3°C - 5°C lower than the primer T<sub>m</sub>).

## Tips for Primer Design

1. Choose C or G as the last base of the 3'-end of the primer.
2. Avoid continuous mismatch at the last 8 bases of the 3'-end of the primer.
3. Avoid hairpin structure at the 3'-end of the primer.
4. The T<sub>m</sub> values of the forward primer and the reverse primer are preferably not more than 1°C, and the T<sub>m</sub> value is preferably adjusted to 55°C-65°C. (Primer Premier 5 is recommended to calculate T<sub>m</sub> value).
5. Primer additional sequence, i.e., unpaired sequence with template, should not be calculated for T<sub>m</sub> value.
6. GC content of the primers should be between 40% - 60%.
7. The overall distribution of primers A, G, C, and T should be as uniform as possible, avoiding the use of high GC/AT content region.
8. Avoid a complementary sequence of more than 5 bases inside of the primer or between the two primers. Avoid a complementary sequence of more than 3 bases at the 3'-end of the two primers.
9. Use the [NCBI BLAST](#) function to search for primer specificity to avoid non-specific amplification.



ISO 9001: 2015



Vazyme Biotech Co., Ltd  
www.vazyme.com

Order: [global@vazyme.com](mailto:global@vazyme.com)

Support: [global@vazyme.com](mailto:global@vazyme.com)

**For research use only, not for use in diagnostic procedures.**